ALBUMIN

Cat. No.	Pack Name	Packaging (Content)
BLT00001	ALB 250	R1: 5 × 50 mL, R2 standard 1 × 3 mL, instruction for use

INTENDED USE

The kit is intended for *in vitro* photometric quantitative determination of albumin in human serum and plasma on various automatic systems. Intended for screening and monitoring of the level of albumin and diagnosis of hypoalbuminemia or hyperalbuminemia. For professional use in clinical laboratories only.

CLINICAL SIGNIFICANCE

Albumin is a carbohydrate-free protein, which constitutes 55–65 % of total plasma protein. It maintains plasma oncotic pressure, and is also involved in the transport and storage of a wide variety of ligands and is a source of endogenous amino acids. Albumin binds and solubilizes various compounds, e.g. bilirubin, calcium and long-chain fatty acids. Furthermore, albumin is capable of binding toxic heavy metal ions as well as numerous pharmaceuticals, which is the reason why lower albumin concentrations in blood have a significant effect on pharmacokinetics.

Increased levels are observed in

- Dehydration due to reduced plasma water content.
- Stasis during venipuncture which causes fluid to escape into the extravascular compartment.

Decreased levels are observed in

- Excessive protein loss (mainly albumin) - from kidney, skin or intestine.
- Decreased synthesis due to dietary, hepatic disease or malabsorption.
- Increased catabolism in fever, untreated diabetes mellitus and hypertension.

PRINCIPLE

At a pH value of 4.2, albumin displays a sufficiently cationic character to be able to bind with bromocresol green (BCG), an anionic dye, to form a blue-green complex:

Albumin + BCG $\longrightarrow$ albumin-BCG complex

Absorbance of the blue-green complex measured at 540–630 nm is proportional to the concentration of albumin in the sample.

REAGENT DESCRIPTION AND COMPOSITION

R1		R2 standard	
Succinate buffer (pH 4.2)	100 mmol/L	Albumin	See bottle label
Bromocresol green	0.21 mmol/L		
Sodium azide	0.5 g/L		

COMPOSITION OF REACTION MIXTURE

Succinate buffer (pH 4.2)	99 mmol/L
Bromocresol green	0.208 mmol/L
Sodium azide	0.495 g/L

REAGENT PREPARATION

Reagents are liquid, ready to use.

MATERIAL REQUIRED BUT NOT PROVIDED WITH THE DEVICE

Any instrument with temperature control of 37 ±0.5 °C that is capable of reading absorbance at 540–630 nm may be used, general laboratory equipment.

XL MULTICAL 4×3, Cat. No. XSYS0034
XL MULTICAL 10×3, Cat. No. XSYS0122

ERBA NORM 4×5, Cat. No. BLT00080
ERBA NORM 10×5, Cat. No. XSYS0123

ERBA PATH 4×5, Cat. No. BLT00081
ERBA PATH 10×5, Cat. No. XSYS0124

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C. Reagents are ready to use. After opening, reagent is stable until expiry date at 2–8 °C if stored at appropriate conditions, closed carefully, protected from light and without any contamination.

SPECIMEN COLLECTION AND HANDLING

It is recommended to follow ISO 15189 and laboratory instruction.

For specimen collection and preparation only use suitable tubes or collection containers.

Only the specimens listed below were tested and found acceptable.

Serum.

Plasma: Li-heparin and K₂-EDTA plasma.

Do not use fluoride plasma.

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer. Centrifuge samples containing precipitates before performing the assay.

See the Limitations and Interferences section for details about possible sample interferences.

Stability in serum / plasma [®] :	2.5 months at	20–25 °C
	5 months at	4–8 °C
	4 months at	-20 °C

Discard contaminated specimens.

CALIBRATION

Calibration with calibrator XL MULTICAL is recommended.

2 point calibration (blank and calibrator); distilled water is recommended as blank

Calibration frequency: it is recommended to do a calibration

- after reagent lot change
- as required by internal quality control procedures

QUALITY CONTROL

For quality control ERBA NORM and ERBA PATH are recommended.

The control intervals and limits should be adapted according to each individual laboratory's requirements. Values obtained should fall within the defined intervals. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

TRACEABILITY

This method, calibrator XL MULTICAL, R2 standard and controls ERBA NORM and ERBA PATH have been standardized against the reference material BCR470/CRM470[®].

ASSAY PROCEDURE

Wavelength: 578 (540–630) nm

Cuvette: 1 cm

	Reagent blank	Calibrator (Standard)	Sample
Reagent 1	1.000 mL	1.000 mL	1.000 mL
Sample	–	–	0.010 mL
Calibrator	–	0.010 mL	–
Distilled water	0.010 mL	–	–

Mix and incubate 1–5 min. at 37 °C. Measure absorbance of the sample A_{sam} and calibrator (standard) A_{cal} against reagent blank.

CALCULATION

$$\text{Albumin (g/dL)}^* = \frac{A_{\text{sam}}}{A_{\text{cal}}} \times C_{\text{cal}} \quad C_{\text{cal}} = \text{calibrator (standard) concentration}$$

ASSAY PARAMETERS FOR PHOTOMETERS

Mode	End Point	Normal High (g/dL)	5.2
Wavelength 1 (nm)	600	Linearity Low (mg/dL)	0.07
Sample Volume (μL)	5/10	Linearity High (mg/dL)	6.50
Working Reagent Volume (μL)	500/1000	Concentration of Standard	See bottle label
Incubation time (min.)	1	Blank with	Reagent
Reaction temperature (°C)	37	Absorbance limit (max.)	0.4
Reaction direction	Increasing	Units	g/dL
Normal Low (g/dL)	3.5		

UNIT CONVERSION

g/dL × 10 = g/L

EXPECTED VALUES[®]

In Serum:

Adults	20–60 years:	3.5–5.2 g/dL	Children	0–4 days	2.8–4.4 g/dL
	60–90 years:	3.2–4.6 g/dL		4 days–14 years	3.8–5.4 g/dL
				14–18 years	3.2–4.5 g/dL

It is recommended that each laboratory verifies this range or derives reference interval for the population it serves.

ANALYTICAL PERFORMANCE

Data contained within this section is representative for performance on ERBA XL-640 automatic system. Data obtained in your laboratory may differ from these values.

Limit of quantification: 0.07 g/dL

Limit of quantification represents the lowest measurable analyte level. It is calculated as the determined activity of diluted sample to have CV <20 % (n = 30).

Linearity: 6.50 g/dL

Linearity is the highest measured activity with recovery within ±10 % from theoretical value.

Precision:

Precision was determined by using controls in an internal protocol with repeatability (n = 20) and intermediate precision (2 aliquots per run, 2 run per day, 20 days). The following results were obtained:

Repeatability	Mean (g/dL)	SD (g/dL)	CV (%)
Sample 1	3.46	0.016	0.46
Sample 2	4.85	0.025	0.52

Intermediate precision	Mean (g/dL)	SD (g/dL)	CV (%)
Sample 1	3.43	0.079	2.31
Sample 2	4.86	0.066	1.35

Accuracy

Two different validated control materials were used. Determined bias is and -6.1 % at the target value 4.48 g/dL and -8.6 % at the target value 5.26 g/dL.

Comparison

A comparison between XL-640 automatic system ALBUMIN (y) and a commercially available test (x) using 146 samples gave following results:

Linear regression:

$$y = 0.940x + 0.132 \text{ g/dL} \quad r = 0.977$$

Passing-Bablok[®]:

$$y = 0.962x + 0.064 \text{ g/dL} \quad r = 0.963$$

Interferences

Criterion: Recovery within ±10 % of initial value of albumin concentration in the sample (serum) without interfering substance.

Following substances do not interfere: haemoglobin up to 4.5 g/L, bilirubin up to 40 mg/dL, triglycerides up to 650 mg/dL.

Limitations:

- Deteriorated reagents (e.g. exceeding the storage temperature) may give incorrect results. Maximum allowable absorbance of the reagent blank measured at 600 nm against the distilled water is 0.4.
- High concentration of haemoglobin, bilirubin and triglycerides in sample can interfere with determination of albumin. See paragraph Interferences.

Drugs: No interference was found at therapeutic concentrations using common drug panels[®].

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

Hazards identification in accordance with Regulation (EC) No 1272/2008

R1

UF1: Y8RN-RJU5-0E7S-53CM

Warning

Hazard statement:

H319 Causes serious eye irritation.

Precautionary statement:

P280 Wear protective gloves/protective clothing/eye protection.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P337 + P313 If eye irritation persists: Get medical advice/attention.

R2

Reagent is not classified as dangerous.

WASTE MANAGEMENT

Please refer to local legal requirements.

Erba Lachema s.r.o., Karásek 2219/1d, 621 00 Brno, CZ
e-mail: diagnostics@erba.com, www.erba.com

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1. Dumas BT, Watson WA, Biggs HG. Albumin standards and the measurement of serum albumin with bromocresol green. Clin Chim Acta 1971;31:87-96.
2. Leonard, P. L., Persaud, J., Motwani, R.: Clin. Chim. Acta 35, 409, 1971.
3. Grant GH, Silverman LM, Christenson RH. Amino acids and proteins. In: Tietz NW, ed. Fundamentals of Clinical Chemistry, 3rd edition Philadelphia, PA: WB Saunders 1987:328-330.
4. Marshall WJ, ed. Illustrated Textbook of Clinical Chemistry, 3rd ed. London: Gower Medical Publishing 1989;207-218.
5. Use of Anticoagulants in Diagnostic Laboratory Investigations. WHO Publication WHO/DIL/LAB/99.1 Rev. 2: Jan 2002.
6. Baudner S, Bienvenu J, Blirup-Jensen S, et al. The certification of a matrix reference material for immunochemical measurement of 14 human serum proteins CRM470. Report EUR 15243 EN 1993;1-186.
7. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. Burtis, C.A., Ashwood, E.R., Bruns, D.E.; 5th edition, WB Saunders Comp., 2012.
8. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11): 783-790.
9. Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 2001;38:376-385

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