ALKALINE PHOSPHATASE

Cat. No.	Pack Name	Packaging (Content)
BLT00003	ALP AMP 150	R1: 4 × 30 mL, R2: 1 × 30 mL instruction for use
BLT00004	ALP AMP 500	R1: 4 × 100 mL, R2: 1 × 100 mL instruction for use

INTENDED USE


The kit is intended for *in vitro* photometric quantitative determination of alkaline phosphatase in human serum and plasma on various automatic systems. In combination with other parameters it is intended for screening, monitoring and diagnosis of liver disease and bone disease. For professional use in clinical laboratories only.

CLINICAL SIGNIFICANCE

Human ALP consists of a group of enzymes which hydrolyse phosphates at an alkaline pH. Alkaline phosphatase in serum consists of four structural genotypes: the liver-bone-kidney type, the intestinal type, the placental type and the variant from the germ cells. It occurs in osteoblasts, hepatocytes, leukocytes, the kidneys, spleen, placenta, prostate and the small intestine. The liver-bone-kidney type is particularly important. A rise in the alkaline phosphatase occurs with all forms of cholestasis, particularly with obstructive jaundice. It is also elevated in diseases of the skeletal system, such as Paget's disease, hyperparathyroidism, rickets and osteomalacia or fractures. A considerable rise in the alkaline phosphatase activity is sometimes seen in children and juveniles. It is caused by increased osteoblast activity following accelerated bone growth.

PRINCIPLE

The products are based on the principle of photometric determination according to IFCC recommendation¹. This method utilises 4-nitrophenyl phosphate as the substrate. Under optimised conditions ALP present in the sample catalyses the following reaction:

At the pH of the reaction, 4-nitrophenol has an intense yellow colour. The reagent also contains a metal ion buffer system to ensure that optimal concentrations of Zinc and Magnesium are maintained. The metal ion buffer can also chelate other potentially inhibitory ions which may be present. The reaction is monitored by measuring the rate of increase in absorbance at 405–430 nm which is proportional to the activity of ALP in the serum.

REAGENT DESCRIPTION AND COMPOSITION

R1		R2	
AMP buffer, pH 10.5	434 mmol/L	AMP buffer, pH 9.5	50 mmol/L
Magnesium acetate	2.48 mmol/L	Magnesium acetate	0.28 mmol/L
Zinc sulfate	1.24 mmol/L	Zinc sulfate	0.14 mmol/L
HEDTA	2.48 mmol/L	HEDTA	0.28 mmol/L
Sodium azide	0.9 g/L	p-nitrophenyl phosphate	81.6 mmol/L
		Sodium azide	0.9 g/L

COMPOSITION OF REACTION MIXTURE

AMP buffer	350 mmol/L
Magnesium acetate	2.00 mmol/L
Zinc sulfate	1.00 mmol/L
HEDTA	2.00 mmol/L
p-nitrophenyl phosphate	16.0 mmol/L
Sodium azide	0.9 g/L

REAGENT PREPARATION

Reagents are liquid, ready to use. For monoreagent method, prepare working reagent by mixing of 4 portions of reagent R1 with 1 portion of reagent R2.

MATERIAL REQUIRED BUT IS NOT PROVIDED WITH THE DEVICE

Any instrument with temperature control of 37 ±0.5 °C that is capable of reading absorbance at 405–430 nm may be used, general laboratory equipment.
 XL MULTICAL 4×3, Cat. No. XSYS0034
 XL MULTICAL 10×3, Cat. No. XSYS0122
 ERBA NORM 4×5, Cat. No. BLT00080
 ERBA NORM 10×5, Cat. No. XSYS0123
 ERBA PATH 4×5, Cat. No. BLT00081
 ERBA PATH 10×5, Cat. No. XSYS0124

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.

Two reagents method – substrate start

Reagents are ready to use. After opening, reagents are stable until expiry date at 2–8 °C if stored at appropriate conditions, closed carefully, protected from light and without any contamination.

Monoreagent method – sample start

Stability of working reagent: 1 week at 15–25 °C in dark
 4 weeks at 2–8 °C in dark

SPECIMEN COLLECTION AND HANDLING

It is recommended to follow ISO 15189 and laboratory instruction.
 For specimen collection and preparation only use suitable tubes or collection containers.
 Only the specimens listed below were tested and found acceptable.

Serum.

Plasma: Li-heparin plasma.

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer. Centrifuge samples containing precipitates before performing the assay.
 See the limitations and interferences section for details about possible sample interferences.

Stability in serum / plasma²:	7 days at	20–25 °C
	7 days at	4–8 °C
	2 months at	–20 °C

Discard contaminated specimens.

CALIBRATION

Calibration with calibrator XL MULTICAL is recommended.
 2 point calibration (blank and calibrator); distilled water is recommended as blank
 Calibration frequency: it is recommended to do a calibration
 • after reagent lot change
 • as required by internal quality control procedures

QUALITY CONTROL

For quality control ERBA NORM and ERBA PATH are recommended. The control intervals and limits should be adapted according to each individual laboratory's requirements. Values obtained should fall within the defined intervals. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

TRACEABILITY

This method, calibrator XL MULTICAL and controls ERBA NORM and ERBA PATH have been standardized to original method IFCC¹.

ASSAY PROCEDURE

Wavelength: 420 (405–430) nm

Cuvette: 1 cm

Two reagents method – substrate start

	Reagent blank	Calibrator	Sample
Reagent 1	0.800 mL	0.800 mL	0.800 mL
Sample	–	–	0.020 mL
Calibrator	–	0.020 mL	–
Distilled water	0.020 mL	–	–

Mix and after 5 min. incubation (at 37 °C) add:

Reagent 2	0.200 mL	0.200 mL	0.200 mL
-----------	----------	----------	----------

Mix, incubate 1 min. at 37 °C and then measure the initial absorbance of calibrator and sample against reagent blank. Measure the absorbance change exactly after 1, 2 and 3 min. Calculate 1 minute absorbance change (ΔA/min).

Monoreagent method – sample start

	Reagent blank	Calibrator	Sample
Working reagent	1.000 mL	1.000 mL	1.000 mL
Sample	–	–	0.020 mL
Calibrator	–	0.020 mL	–
Distilled water	0.020 mL	–	–

Mix, incubate 1 min. at 37 °C and then measure the initial absorbance of calibrator and sample against reagent blank. Measure the absorbance change exactly after 1, 2 and 3 min. Calculate 1 minute absorbance change (ΔA/min).

CALCULATION

- ALP (U/L) = $\frac{\Delta A_{\text{sam}} / \text{min.}}{\Delta A_{\text{cal}} / \text{min.}} \times C_{\text{cal}}$ C_{cal} = calibrator concentration
- Using factor:
 ALP (U/L) = f × ΔA/min
 f = 2764 (at 405 nm) f = factor

ASSAY PARAMETERS FOR PHOTOMETERS

Mode	Kinetic	Reaction direction	Increasing
Wavelength (nm)	405	Normal Low U/L	42
Sample Volume (μL)	10/20	Normal High U/L	128
Working Reagent Volume (μL)	500/1000	Linearity Low U/L	4.94
Lag time (sec.)	60	Linearity High U/L	1617
Kinetic interval (sec.)	60	Blank with	Reagent
No. of readings	3	Absorbance limit (max.)	1.4
Kinetic factor	2764	Units	U/L
Reaction temperature (°C)	37		

UNIT CONVERSION

U/L × 0.0167 = μkat/L

EXPECTED VALUES³

At 37 °C

Females: 4–15 years:	54–369 U/L	Males: 1–12 years:	54–369 U/L
20–50 years:	42–98 U/L	20–50 years:	53–128 U/L
≥60 years:	53–141 U/L	≥60 years:	56–119 U/L

It is recommended that each laboratory verifies this range or derives reference interval for the population it serves.

ANALYTICAL PERFORMANCE

Data contained within this section are representative for performance on ERBA XL-640 automatic system. Data obtained in your laboratory may differ from these values.

Limit of quantification:

4.94 U/L
 Limit of quantification represents the lowest measurable analyte level. It is calculated as the determined activity of diluted sample to have CV <20 % (n = 30).

Linearity:

1617 U/L
 Linearity is the highest measured activity with recovery within ±10 % from theoretical value.

Precision:

Precision was determined by using controls in an internal protocol with repeatability (n = 20) and intermediate precision (2 aliquots per run, 2 run per day, 20 days). The following results were obtained:

Repeatability	Mean (U/L)	SD (U/L)	CV (%)
Sample 1	125	0.55	0.44
Sample 2	374	1.41	0.38

Intermediate precision	Mean (U/L)	SD (U/L)	CV (%)
Sample 1	130	3.77	2.90
Sample 2	392	7.86	2.01

Accuracy

Two different validated control materials were used. Determined bias is –4.1 % at the target value 163.8 U/L and –1.6 % at the target value 392.1 U/L.

Comparison

A comparison between XL-640 automatic system ALP (y) and a commercially available test (x) using 47 samples gave following results:

Linear regression:

$$y = 1.049x - 0.18 \text{ U/L} \quad r = 1.000$$

Passing-Bablok¹:

$$y = 1.048x - 0.03 \text{ U/L} \quad r = 1.000$$

Interferences

Criterion: Recovery within ±10 % of initial value of an alkaline phosphatase activity in the sample without interfering substance.

Following substances do not interfere: haemoglobin up to 9 g/L, bilirubin up to 40 mg/dL, triglycerides up to 850 mg/dL.

Limitations:

- Deteriorated reagents (e.g. exceeding the storage temperature) may give incorrect results. Maximum allowable absorbance of the reagent blank measured at 420 nm against the distilled water is 1.0.
- High concentration of haemoglobin, bilirubin and triglycerides in sample can interfere with determination of ALP. See paragraph Interferences.

Drugs: No interference was found at therapeutic concentrations using common drug panels⁸.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

Hazards identification in accordance with Regulation (EC) No 1272/2008

R1

The reagent is not classified as dangerous.

R2

UFI: 5FAQ-SJAM-HE7A-6Q8K

Warning

Hazard statement:

- H315 Causes skin irritation.
- H319 Causes serious eye irritation.

Precautionary statement:

- P280 Wear protective gloves / protective clothing / eye protection.
- P302 + P352 IF ON SKIN: Wash with plenty of water.
- P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

WASTE MANAGEMENT

Please refer to local legal requirements.

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Table with 4 columns: Reagent, Blank reagent, Calibrator, Sample. Rows show the procedure for the two-reagent method.

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Table with 4 columns: Working reagent, Blank reagent, Calibrator, Sample. Rows show the procedure for the mono-reagent method.

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Table with 4 columns: Mode, Kinetic, Reaction direction, Amplification. Rows list various parameters for the photometer analysis.

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Table showing expected values for males and females across different age groups (4-15, 20-50, ≥60 years).

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Table with 4 columns: Repeatability, Mean (Ед/л), SD (Ед/л), CV (%). Rows show repeatability and intermediate precision results for two samples.

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For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.
Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

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1. IFCC method for the measurement of ALP J. Clin. Chem. Clin. Biochem. 1983; 21: 731-48.
2. Guder WG, Narayanan S, Wisser H, et al. List of Analytes; Preanalytical Variables. Brochure in: Samples: From the Patient to the Laboratory. Darmstadt: GIT-Verlag 1996.
3. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. Burtis, C.A., Ashwood, E.R., Bruns, D.E.; 5th edition, WB Saunders Comp., 2012.
4. Zilva JF, Pannall PR, "Plasma Enzymes in Diagnosis" in Clinical Chemistry in Diagnosis and Treatment. Lloyd London 1979: Chapter 15: 343.
5. Young DS. Effects of Drugs on Clinical Laboratory Tests. Third Edition 1990: 3: 19-25.
6. Kaplan and Pesce (Eds.) Clinical Chemistry, Theory analysis and correlation. Second Edition. CV Mosby Co. 1989.
7. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11):783-790.
8. Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 2001;38:376-385

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