AST/GOT

Cat. No.	Pack Name	Packaging (Content)
BLT00050	AST/GOT 250	R1: 4 × 50 mL, R2: 1 × 50 mL, instruction for use
BLT00051	AST/GOT 500	R1: 4 × 100 mL, R2: 1 × 100 mL, instruction for use

INTENDED USE


The kit is intended for *in vitro* photometric quantitative determination of aspartate aminotransferase in human serum and plasma on various automatic systems. In combination with other parameters it is intended for screening, monitoring and diagnosis of liver diseases. For professional use in clinical laboratories only.

CLINICAL SIGNIFICANCE

Aspartate aminotransferase (AST/GOT) is widely distributed with high concentrations in the heart, liver, skeletal muscle, kidney and erythrocytes. Damage or disease to any of these tissues such as myocardial infarction, viral hepatitis, liver necrosis, cirrhosis and muscular dystrophy may result in raised levels of AST/GOT. Following myocardial infarction, serum AST is elevated and reaches a peak two days after onset. In patients undergoing renal dialysis or those with vitamin B6 deficiency, serum AST/GOT may be decreased. The apparent reduction in AST may be related to decreased pyridoxal phosphate, the prosthetic group for AST/GOT, resulting in an increase in the ratio of apoenzyme to holoenzyme. Addition of pyridoxal-5-phosphate, recommended by IFCC¹, stabilizes the activity of AST/GOT and avoids falsely low values in samples containing insufficient endogenous pyridoxal-5-phosphate, e.g. from patients with myocardial infarction, liver disease and intensive care patients.

PRINCIPLE

This assay follows the recommendations of the IFCC¹ without pyridoxal phosphate. AST in the sample catalyzes the transfer of an amino group between L-aspartate and 2-oxoglutarate to form oxaloacetate and L-glutamate. The oxaloacetate then reacts with NADH, in the presence of malate dehydrogenase (MDH), to form NAD⁺. The series of reactions involved in the assay system is as follows:

The reaction is monitored by measuring the rate of decrease in absorbance at 340 nm due to the oxidation of NADH. The rate of oxidation of NADH is proportional to the activity of AST/GOT present in the sample.

Addition of lactate dehydrogenase to the reagent is necessary to achieve rapid and complete reduction of endogenous pyruvate so that it does not interfere with the assay.

REAGENT DESCRIPTION AND COMPOSITION

R1		R2	
Tris buffer (pH 7.8)	110 mmol/L	CAPSO	20 mmol/L
L-aspartic acid	340 mmol/L	2-oxoglutarate	85 mmol/L
LDH	≥2250 U/L	NADH	1.05 mmol/L
MDH	≥810 U/L	Sodium azide	0.95 g/L
Sodium chloride	200 mmol/L		
Sodium azide	0.95 g/L		

COMPOSITION OF REACTION MIXTURE

Tris buffer	81.5 mmol/L
L-aspartic acid	252 mmol/L
LDH	≥1667 U/L
MDH	≥600 U/L
CAPSO	3.7 mmol/L
2-oxoglutarate	15.7 mmol/L
NADH	0.19 mmol/L
Sodium chloride	148 mmol/L
Sodium azide	0.88 g/L

REAGENT PREPARATION

Reagents are liquid, ready to use. For monoreagent method, prepare working reagent by mixing of 4 portions of reagent R1 with 1 portion of reagent R2.

MATERIAL REQUIRED BUT NOT PROVIDED WITH THE DEVICE

Any instrument with temperature control of 37 ±0.5 °C that is capable of reading absorbance at 340 nm may be used, pyridoxal-5-phosphate (PDP 50, Cat. No. 10003294), general laboratory equipment.
XL MULTICAL 4×3, Cat. No. XSYS0034
XL MULTICAL 10×3, Cat. No. XSYS0122
ERBA NORM 4×5, Cat. No. BLT00080
ERBA NORM 10×5, Cat. No. XSYS0123
ERBA PATH 4×5, Cat. No. BLT00081
ERBA PATH 10×5, Cat. No. XSYS0124

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.

Two reagents method – substrate start

Reagents are ready to use. After opening, reagents are stable until expiry date at 2–8 °C if stored at appropriate conditions, closed carefully, protected from light and without any contamination.

Monoreagent method – sample start

Stability of working reagent: 5 days at 15–25 °C
3 weeks at 2–8 °C

SPECIMEN COLLECTION AND HANDLING

It is recommended to follow ISO 15189 and laboratory instruction.

For specimen collection and preparation only use suitable tubes or collection containers.

Only the specimens listed below were tested and found acceptable.

Serum.

Plasma: Li-heparin and K₂-EDTA plasma.

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer. Centrifuge samples containing precipitates before performing the assay.

See the limitations and interferences section for details about possible sample interferences.

Stability in serum / plasma:	4 days at 15–25 °C
	7 days at 2–8 °C
	3 months at -20 °C

Discard contaminated specimens.

CALIBRATION

Calibration with calibrator XL MULTICAL is recommended.
2 point calibration (blank and calibrator); distilled water is recommended as blank
Calibration frequency: it is recommended to do a calibration
• after reagent lot change
• as required by internal quality control procedures

QUALITY CONTROL

For quality control ERBA NORM and ERBA PATH are recommended.

The control intervals and limits should be adapted according to each individual laboratory's requirements. Values obtained should fall within the defined intervals. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

TRACEABILITY

This method, calibrator XL MULTICAL and controls ERBA NORM and ERBA PATH have been standardized to IFCC recommended method¹.

ASSAY PROCEDURE

Wavelength: 340 nm
Cuvette: 1 cm

Two reagents method – substrate start

	Reagent blank	Calibrator	Sample
Reagent 1	0.800 mL	0.800 mL	0.800 mL
Sample	–	–	0.100 mL
Calibrator	–	0.100 mL	–
Distilled water	0.100 mL	–	–

Mix and after 5 min. incubation (at 37 °C) add:

Reagent 2	0.200 mL	0.200 mL	0.200 mL
-----------	----------	----------	----------

Mix, incubate 1 min. at 37 °C and then measure the initial absorbance of calibrator and sample against reagent blank. Measure the absorbance change exactly after 1, 2 and 3 min. Calculate 1 minute absorbance change (ΔA/min).

Monoreagent method – sample start

	Reagent blank	Calibrator	Sample
Working reagent	1.000 mL	1.000 mL	1.000 mL
Sample	–	–	0.100 mL
Calibrator	–	0.100 mL	–
Distilled water	0.100 mL	–	–

Mix, incubate 1 min. at 37 °C and then measure the initial absorbance of calibrator and sample against reagent blank. Measure the absorbance change exactly after 1, 2 and 3 min. Calculate 1 minute absorbance change (ΔA/min).

CALCULATION

$$1. \text{AST/GOT (U/L)} = \frac{\Delta A_{\text{sam}}/\text{min.}}{\Delta A_{\text{cal}}/\text{min.}} \times C_{\text{cal}} \quad C_{\text{cal}} = \text{calibrator concentration}$$

$$2. \text{Using factor:} \quad \text{AST/GOT} = f \times \Delta A/\text{min} \quad f = \text{factor}$$

Substrate start: Factor at 340 nm 2143

Sample start: Factor at 340 nm 1745

ASSAY PARAMETERS FOR PHOTOMETERS

Mode	Kinetic	Reaction direction	Decreasing
Wavelength (nm)	340	Normal Low U/L	0
Sample Volume (μL)	50/100	Normal High U/L	35
Working Reagent Volume (μL)	500/1000	Linearity Low U/L	2.13
Lag time (sec.)	60	Linearity High U/L	418
Kinetic interval (sec.)	60	Blank with	Water
No. of readings	3	Absorbance limit (min.)	1.1
Kinetic factor	1745	Units	U/L
Reaction temperature (°C)	37		

UNIT CONVERSION

U/L × 0.0167 = μkat/L

EXPECTED VALUES⁵

At 37 °C

Females <31 U/L

Males <35 U/L

It is recommended that each laboratory verifies this range or derives reference interval for the population it serves.

ANALYTICAL PERFORMANCE

Data contained within this section is representative for performance on ERBA XL-640 automatic system. Data obtained in your laboratory may differ from these values.

Limit of quantification: 2.13 U/L

Limit of quantification represents the lowest measurable analyte level. It is calculated as the determined activity of diluted sample to have CV <20 % (n = 30).

Linearity: 418 U/L
Linearity is the highest measured activity with recovery within ±10 % from theoretical value.

Precision:

Precision was determined by using controls in an internal protocol with repeatability (n = 20) and intermediate precision (2 aliquots per run, 2 run per day, 20 days). The following results were obtained:

Repeatability	Mean (U/L)	SD (U/L)	CV (%)	Intermediate precision	Mean (U/L)	SD (U/L)	CV (%)
Sample 1	42.4	0.60	1.41	Sample 1	42.5	1.15	2.71
Sample 2	153.3	1.13	0.74	Sample 2	151.7	2.06	1.36

Accuracy

Two different validated control materials were used. Determined bias is -4.8 % at the target value 146.6 U/L and -4.2 % at the target value 197.2 U/L.

Comparison

A comparison between XL-640 automatic system AST/GOT with PDP (y) and a commercially available test (x) using 48 samples gave following results:

Linear regression: y = 1.025x - 2.74 U/L r = 0.999

Passing-Bablok⁶: y = 1.000x - 1.80 U/L r = 0.977

Interferences

Criterion: Recovery within ±10 % of initial value of AST/GOT activity in the sample without interfering substance.
Following substances do not interfere: bilirubin up to 40 mg/dL, triglycerides up to 850 mg/dL. Haemolysis interferes due to AST/GOT activity from erythrocytes.

Drugs: No interference was found at therapeutic concentrations using common drug panels⁷.

Limitations:

- Deteriorated reagents (e.g. exceeding the storage temperature) may give incorrect results. Minimal allowable absorbance of the reagent blank measured at 340 nm against the distilled water is 1.1.
- High concentration of haemoglobin, bilirubin and triglycerides in sample can interfere with determination of AST/GOT. See paragraph Interferences.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

Hazards identification in accordance with Regulation (EC) No 1272/2008

R1

UFI: YMDF-7EN4-1H7J-XF4E

Warning

Contains: aspartic acid, tris(hydroxymethyl)aminomethane, sodium hydroxide

Hazard statement:

H315 Causes skin irritation.

H319 Causes serious eye irritation.

Precautionary statement:

P280 Wear protective gloves / protective clothing / eye protection.

P302 + P352 IF ON SKIN: Wash with plenty of water and soap.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P337 + P313 If eye irritation persists: Get medical advice/attention.

R2

The reagent is not classified as dangerous.

WASTE MANAGEMENT

Please refer to local legal requirements.

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- Schumann G, Bonora R, Ceriotti F, Féraud G et al. IFCC primary reference procedure for the measurement of catalytic activity concentrations of enzymes at 37 °C. Part 5: Reference procedure for the measurement of catalytic concentration of aspartate aminotransferase. Clin Chem Lab Med 2002;40:725-33.
- Thomas L. Alanine aminotransferase (ALT), Aspartate aminotransferase (AST). In: Thomas L, editor, Clinical Laboratory Diagnostics, 1st ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998, p. 55-65.
- Bergmeyer HU, Horder M, Rej R, Approved Recommendation (1985) on IFCC Methods for the Measurement of Catalytic Concentration of Enzymes, L. Clin. Chem. Clin. Biochem. 1986; 24: 481-495
- Use of Anticoagulants in Diagnostic Laboratory Investigations. WHO Publication WHO/DIL/LAB/99.1 Rev. 2: Jan 2002.
- Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. Burtis, C.A., Ashwood, E.R., Bruns, D.E.; 5th edition, WB Saunders Comp., 2012.
- Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11):783-790.
- Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 2001;38:376-385

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Dátum revízie: 25. 9. 2025

LITERATÚRA

1. Schumann G, Bonora R, Ceriotti F, Férard G et al. IFCC primary reference procedure for the measurement of catalytic activity concentrations of enzymes at 37 °C. Part 5: Reference procedure for the measurement of catalytic concentration of aspartate aminotransferase. Clin Chem Lab Med 2002;40:725-33.
2. Thomas L. Alanine aminotransferase (ALT), Aspartate aminotransferase (AST). In: Thomas L, editor, Clinical Laboratory Diagnostics, 1st ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998, p. 55-65.
3. Bergmeyer HU, Horder M, Rej R, Approved Recommendation (1985) on IFCC Methods for the Measurement of Catalytic Concentration of Enzymes, L. Clin. Chem. Clin. Biochem. 1986; 24: 481-495
4. Use of Anticoagulants in Diagnostic Laboratory Investigations. WHO Publication WHO/DIL/LAB/99.1 Rev. 2: Jan 2002.
5. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. Burtis, C.A., Ashwood, E.R., Bruns, D.E.; 5th edition, WB Saunders Comp., 2012.
6. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11):783-790.
7. Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 2001;38:376-385

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Katalógové číslo

Číslo šarže

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eIFU:
www.erba.com

Diagnostický zdravotnícký prostriedok *in vitro*

Výrobca

Obmedzenie teploty

Obsah

QUALITY SYSTEM CERTIFIED
ISO 13485, IVDR

Erba Lachema s.r.o., Karásek 2219/1d, 621 00 Brno, CZ
e-mail: diagnostics@erba.com, www.erba.com
CC/IFU/007/25/A

Dátum revízie: 25. 9. 2025