ALT/GPT

Cat. No.	Pack Name	Packaging (Content)
BLT00052	ALT/GPT 250	R1: 4 × 50 mL, R2: 1 × 50 mL, instruction for use
BLT00053	ALT/GPT 500	R1: 4 × 100 mL, R2: 1 × 100 mL, instruction for use

INTENDED USE

The kit is intended for *in vitro* photometric quantitative determination of alanine aminotransferase in human serum and plasma on various automatic systems. In combination with other parameters it is intended for screening, monitoring and diagnosis of liver diseases. For professional use in clinical laboratories only.

CLINICAL SIGNIFICANCE

Alanine aminotransferase (ALT/GPT) is present in high concentration in liver and to a lesser extent in kidney, heart, skeletal muscle, pancreas, spleen and lung. Increased levels of ALT/GPT however is generally a result of liver disease associated with some degree of hepatic necrosis such as cirrhosis, viral or toxic hepatitis and obstructive jaundice.


Characteristically ALT/GPT is generally higher than aspartate aminotransferase (AST/GOT) in acute viral or toxic hepatitis, whereas for most patients with chronic hepatic disease, ALT/GPT levels are generally lower than AST/GOT levels. Elevated ALT/GPT levels have also been found in extensive trauma and muscle disease, circulatory failure with shock, hypoxia, myocardial infarction and haemolytic disease.

Although both sera AST/GOT and ALT/GPT become elevated whenever disease processes affect liver cell integrity, ALT/GPT is more liver-specific enzyme. Moreover, elevations of ALT/GPT activity persist longer than elevations of AST/GOT activity.

In patients with vitamin B6 deficiency, serum aminotransferase activity may be decreased. The apparent reduction in aminotransferase activity may be related to decreased pyridoxal-5-phosphate, the prosthetic group for aminotransferases, resulting in an increase in the ratio of apoenzyme to holoenzyme. Addition of pyridoxal-5-phosphate, recommended by IFCC¹, stabilizes the activity of ALT/GPT and avoids falsely low values in samples containing insufficient endogenous pyridoxal-5-phosphate, e.g. from patients with myocardial infarction, liver disease and intensive care patients.

PRINCIPLE

This assay follows the recommendations of the IFCC¹ without pyridoxal phosphate. ALT/GPT catalyzes the reaction between L-alanine and 2-oxoglutarate. The pyruvate formed is reduced by NADH in a reaction catalyzed by lactate dehydrogenase (LDH) to form L-lactate and NAD⁺. The series of reactions involved in the assay system is as follows:

The reaction is monitored by measuring the rate of decrease in absorbance at 340 nm due the oxidation of NADH. The rate of oxidation of NADH is proportional to the activity of ALT/GPT present in the sample. Endogenous sample pyruvate is rapidly and completely reduced by LDH during initial incubation period to avoid interference during the assay.

REAGENT DESCRIPTION AND COMPOSITION

R1		R2	
Tris buffer (pH 7.5)	137.5 mmol/L	CAPSO	20 mmol/L
L-alanine	709 mmol/L	2-oxoglutarate	85 mmol/L
LDH	≥2000 U/L	NADH	1.05 mmol/L
Sodium chloride	200 mmol/L	Sodium azide	0.95 g/L
Sodium azide	0.95 g/L		

COMPOSITION OF REACTION MIXTURE

Tris buffer	102 mmol/L
L-alanine	525 mmol/L
LDH	≥1480 U/L
CAPSO	3.7 mmol/L
2-oxoglutarate	15.7 mmol/L
NADH	0.19 mmol/L
Sodium chloride	148 mmol/L
Sodium azide	0.88 g/L

REAGENT PREPARATION

Reagents are liquid, ready to use. For monoreagent method, prepare working reagent by mixing of 4 portions of reagent R1 with 1 portion of reagent R2.

MATERIAL REQUIRED BUT NOT PROVIDED WITH THE DEVICE

Any instrument with temperature control of 37 ±0.5 °C that is capable of reading absorbance at 340 nm may be used, pyridoxal-5-phosphate (PDP 50, Cat. No. 10003294), general laboratory equipment.

XL MULTICAL 4×3, Cat. No. XSYS0034
XL MULTICAL 10×3, Cat. No. XSYS0122
ERBA NORM 4×5, Cat. No. BLT00080
ERBA NORM 10×5, Cat. No. XSYS0123
ERBA PATH 4×5, Cat. No. BLT00081
ERBA PATH 10×5, Cat. No. XSYS0124

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.

Two reagents method – substrate start

Reagents are ready to use. After opening, reagents are stable until expiry date at 2–8 °C if stored at appropriate conditions, closed carefully, protected from light and without any contamination.

Monoreagent method – sample start

Stability of working reagent: 5 days at 15–25 °C
4 weeks at 2–8 °C

SPECIMEN COLLECTION AND HANDLING

It is recommended to follow ISO 15189 and laboratory instruction.

For specimen collection and preparation only use suitable tubes or collection containers.

Only the specimens listed below were tested and found acceptable.

Serum.

Plasma: Li-heparin and K₂-EDTA plasma.

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer. Centrifuge samples containing precipitates before performing the assay.

See the limitations and interferences section for details about possible sample interferences.

Stability in serum / plasma^{4,5}:	3 days at	15–25 °C
	7 days at	2–8 °C
	7 days at	–20 °C

Discard contaminated specimens.

CALIBRATION

Calibration with calibrator XL MULTICAL is recommended.

2 point calibration (blank and calibrator); distilled water is recommended as blank

Calibration frequency: it is recommended to do a calibration:

- after reagent lot change
- as required by internal quality control procedures

QUALITY CONTROL

For quality control ERBA NORM and ERBA PATH are recommended.

The control intervals and limits should be adapted according to each individual laboratory's requirements. Values obtained should fall within the defined intervals. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

TRACEABILITY

This method, calibrator XL MULTICAL and controls ERBA NORM and ERBA PATH have been standardized to IFCC recommended method¹.

ASSAY PROCEDURE

Wavelength: 340 nm

Cuvette: 1 cm

Two reagents method – substrate start

	Reagent blank	Calibrator	Sample
Reagent 1	0.800 mL	0.800 mL	0.800 mL
Sample	–	–	0.100 mL
Calibrator	–	0.100 mL	–
Distilled water	0.100 mL	–	–

Mix and after 5 min. incubation (at 37 °C) add:

Reagent 2	0.200 mL	0.200 mL	0.200 mL
-----------	----------	----------	----------

Mix, incubate 1 min. at 37 °C and then measure the initial absorbance of calibrator and sample against reagent blank. Measure the absorbance change exactly after 1, 2 and 3 min. Calculate 1 minute absorbance change (ΔA/min).

Monoreagent method – sample start

	Reagent blank	Calibrator	Sample
Working reagent	1.000 mL	1.000 mL	1.000 mL
Sample	–	–	0.100 mL
Calibrator	–	0.100 mL	–
Distilled water	0.100 mL	–	–

Mix, incubate 1 min. at 37 °C and then measure the initial absorbance of calibrator and sample against reagent blank. Measure the absorbance change exactly after 1, 2 and 3 min. Calculate 1 minute absorbance change (ΔA/min).

CALCULATION

$$1. \text{ ALT/GPT (U/L)} = \frac{\Delta A_{\text{sam}}/\text{min.}}{\Delta A_{\text{cal}}/\text{min.}} \times C_{\text{cal}} \quad C_{\text{cal}} = \text{calibrator concentration}$$

$$2. \text{ Using factor: } \text{ALT/GPT (U/L)} = f \times \Delta A/\text{min} \quad f = \text{factor}$$

Substrate start:
factor at 340 nm 2143

Sample start:
factor at 340 nm 1745

ASSAY PARAMETERS FOR PHOTOMETERS

Mode	Kinetic	Reaction direction	Decreasing
Wavelength (nm)	340	Normal Low U/L	0
Sample Volume (μL)	50/100	Normal High U/L	45
Working Reagent Volume (μL)	500/1000	Linearity Low U/L	2.84
Lag time (sec.)	60	Linearity High U/L	494
Kinetic interval (sec.)	60	Blank with	Water
No. of readings	3	Absorbance limit (min.)	1.1
Kinetic factor	1745	Units	U/L
Reaction temperature (°C)	37		

UNIT CONVERSION

U/L × 0.0167 = μkat/L

EXPECTED VALUES*

At 37 °C

Females <34 U/L

Males <45 U/L

It is recommended that each laboratory verifies this range or derives reference interval for the population it serves.

ANALYTICAL PERFORMANCE

Data contained within this section is representative for performance on ERBA XL-640 automatic system. Data obtained in your laboratory may differ from these values.

Limit of quantification:

2.84 U/L
Limit of quantification represents the lowest measurable analyte level. It is calculated as the determined activity of diluted sample to have CV <20 % (n = 30)

Linearity:

494 U/L
Linearity is the highest measured activity with recovery within ±10 % from theoretical value.

Precision:

Precision was determined by using controls in an internal protocol with repeatability (n = 20) and intermediate precision (2 aliquots per run, 2 run per day, 20 days). The following results were obtained:

Repeatability	Mean (U/L)	SD (U/L)	CV (%)
Sample 1	42.4	0.59	1.38
Sample 2	125.1	1.08	0.86

Intermediate precision	Mean (U/L)	SD (U/L)	CV (%)
Sample 1	42.7	1.13	2.65
Sample 2	123.5	2.10	1.70

Accuracy

Two different validated control materials were used. Determined bias is -4.8 % at the target value 73.3 U/L and -7.2 % at the target value 230.8 U/L.

Comparison

A comparison between XL-640 automatic system ALT/GPT (y) with PDP and a commercially available test (x) using 47 samples gave following results:

Linear regression:

$$y = 1.003x + 0.31 \text{ U/L} \quad r = 0.996$$

Passing-Bablok⁷:

$$y = 1.000x + 0.00 \text{ U/L} \quad r = 0.982$$

Interferences

Criterion: Recovery within ±10 % of initial value of ALT/GPT activity in the sample without interfering substance.

Following substances do not interfere: haemoglobin up to 5 g/L, bilirubin up to 40 mg/dL, triglycerides up to 850 mg/dL.

Drugs: No interference was found at therapeutic concentrations using common drug panels⁸.

Limitations:

- Deteriorated reagents (e.g. exceeding the storage temperature) may give incorrect results. Minimal allowable absorbance of the reagent blank measured at 340 nm against the distilled water is 1.1.
- High concentration of haemoglobin, bilirubin and triglycerides in sample can interfere with determination of ALT/GPT. See paragraph Interferences.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

Hazards identification in accordance with Regulation (EC) No 1272/2008

R1, R2

Reagents of the kit are not classified as dangerous.

WASTE MANAGEMENT

Please refer to local legal requirements.

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- Schumann G, Bonora R, Ceriotti F, Féraud G et al. IFCC primary reference procedure for the measurement of catalytic activity concentrations of enzymes at 37 °C. Part 5: Reference procedure for the measurement of catalytic concentration of alanine aminotransferase. Clin Chem Lab Med 2002;40:718-24.
- Thomas L. Alanine aminotransferase (ALT), Aspartate aminotransferase (AST). In: Thomas L, editor, Clinical Laboratory Diagnostics, 1st ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998, p. 55-65.
- Bergmeyer HU, Horder M, Rej R, Approved Recommendation (1985) on IFCC Methods for the Measurement of Catalytic Concentration of Enzymes, L. Clin. Chem. Clin. Biochem. 1986; 24: 481-495
- Heins M, Heil W, Withold W. Storage of Serum or Whole Blood Samples? Effects of Time and Temperature on 22 Serum Analytes. Eur J Clin Chem Clin Biochem 1995;33:231-238.
- Guder WG, Zawta B et al. The Quality of Diagnostic Samples. 1st ed. Darmstadt: GIT Verlag; 2001; p. 14-5.
- Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. Burtis, C.A., Ashwood, E.R., Bruns, D.E.; 5th edition, WB Saunders Comp., 2012.
- Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11):783-790.
- Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 2001;38:376-385

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Dátum revízie: 25. 9. 2025

LITERATÚRA

1. Schumann G, Bonora R, Ceriotti F, Férard G et al. IFCC primary reference procedure for the measurement of catalytic activity concentrations of enzymes at 37 °C. Part 5: Reference procedure for the measurement of catalytic concentration of alanine aminotransferase. Clin Chem Lab Med 2002;40:718-24.
2. Thomas L. Alanine aminotransferase (ALT), Aspartate aminotransferase (AST). In: Thomas L, editor, Clinical Laboratory Diagnostics, 1st ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998, p. 55-65.
3. Bergmeyer HU, Horder M, Rej R, Approved Recommendation (1985) on IFCC Methods for the Measurement of Catalytic Concentration of Enzymes, L. Clin. Chem. Clin. Biochem. 1986; 24: 481-495
4. Heins M, Heil W, Withold W. Storage of Serum or Whole Blood Samples? Effects of Time and Temperature on 22 Serum Analytes. Eur J Clin Chem Clin Biochem 1995;33:231-238.
5. Guder WG, Zawta B et al. The Quality of Diagnostic Samples. 1st ed. Darmstadt: GIT Verlag; 2001; p. 14-5.
6. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. Burtis, C.A., Ashwood, E.R., Bruns, D.E.; 5th edition, WB Saunders Comp., 2012.
7. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11):783-790.
8. Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 2001;38:376-385

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