ALBUMIN

Cat. No.	Pack Name	Packaging (Content)
XSYS0001	ALB 440	R1: 10 × 44 mL, RFID tag, instruction for use

INTENDED USE

The kit is intended for *in vitro* photometric quantitative determination of albumin in human serum and plasma on automatic systems ERBA XL. Intended for screening and monitoring of the level of albumin and diagnosis of hypoalbuminemia or hyperalbuminemia. For professional use in clinical laboratories only.

CLINICAL SIGNIFICANCE

Albumin is a carbohydrate-free protein, which constitutes 55–65 % of total plasma protein. It maintains plasma oncotic pressure, and is also involved in the transport and storage of a wide variety of ligands and is a source of endogenous amino acids. Albumin binds and solubilizes various compounds, e.g. bilirubin, calcium and long-chain fatty acids. Furthermore, albumin is capable of binding toxic heavy metal ions as well as numerous pharmaceuticals, which is the reason why lower albumin concentrations in blood have a significant effect on pharmacokinetics.

Increased levels are observed in

- Dehydration due to reduced plasma water content.
- Stasis during venipuncture which causes fluid to escape into the extravascular compartment.

Decreased levels are observed in

- Excessive protein loss (mainly albumin) - from kidney, skin or intestine.
- Decreased synthesis due to dietary, hepatic disease or malabsorption.
- Increased catabolism in fever, untreated diabetes mellitus and hypertension.

PRINCIPLE

At a pH value of 4.2, albumin displays a sufficiently cationic character to be able to bind with bromocresol green (BCG), an anionic dye, to form a blue-green complex¹.

Albumin + BCG $\longrightarrow$ albumin-BCG complex

Absorbance of the blue-green complex measured at 600 nm is proportional to the concentration of albumin in the sample.

REAGENT DESCRIPTION AND COMPOSITION

R1	
Succinate buffer (pH 4.2)	100 mmol/L
Bromocresol green	0.21 mmol/L
Sodium azide	0.5 g/L

COMPOSITION OF REACTION MIXTURE

Succinate buffer (pH 4.2)	99 mmol/L
Bromocresol green	0.208 mmol/L
Sodium azide	0.495 g/L

REAGENT PREPARATION

Reagents are liquid, ready to use. Load the number of tests from the RFID tag before using a new kit.

MATERIAL REQUIRED BUT NOT PROVIDED WITH THE DEVICE

XL MULTICAL 4×3, Cat. No. XSYS0034
 XL MULTICAL 10×3, Cat. No. XSYS0122
 ERBA NORM 4×5, Cat. No. BLT00080
 ERBA NORM 10×5, Cat. No. XSYS0123
 ERBA PATH 4×5, Cat. No. BLT00081
 ERBA PATH 10×5, Cat. No. XSYS0124
 Erba XL analysers: XL-200, Cat. No. INS00002
 XL-640, Cat. No. INS00008
 XL-1000, Cat. No. INS00010

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.
 On board stability: min. 60 days if refrigerated (2–10 °C) and not contaminated.

SPECIMEN COLLECTION AND HANDLING

It is recommended to follow ISO 15189 and laboratory instruction.
 For specimen collection and preparation only use suitable tubes or collection containers.
 Only the specimens listed below were tested and found acceptable.
 Serum.

Plasma: Li-heparin and K₂-EDTA plasma.

Do not use fluoride plasma.

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay.

See the Limitations and Interferences section for details about possible sample interferences.

Stability in serum / plasma²:	2.5 months at	20–25 °C
	5 months at	4–8 °C
	4 months at	-20 °C

Discard contaminated specimens.

CALIBRATION

Calibration with calibrator XL MULTICAL is recommended.
 2 point calibration (blank and calibrator); distilled water is recommended as blank
 Calibration frequency: 30 days
 Calibration is needed:

- after reagent lot change
- as required by internal quality control procedures
- calibration interval may be extended based on acceptable verification of calibration by the laboratory

QUALITY CONTROL

For quality control ERBA NORM and ERBA PATH are recommended.
 The control intervals and limits should be adapted according to each individual laboratory's requirements. Values obtained should fall within the defined intervals. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

TRACEABILITY

This method, calibrator XL MULTICAL and controls ERBA NORM and ERBA PATH have been standardized against the reference material BCR470/CRM470⁶.

ASSAY PROCEDURE AND CALCULATION

ERBA XL automatic systems calculate the concentration of each sample. For assay parameters see www.erba.com.

Assay parameters for ERBA XL automatic systems

Assay type	1-Point
Curve type	Linear
Wavelength (prim. / sec.)	600 nm
Reading time	90–180 s after adding of R1
Reaction direction	Increase
Unit	g/dL (g/L)
Reagent volumes	
R1	200 µL
Sample	2 µL

Note: reagents and sample volumes can be different for individual ERBA XL automatic systems depending on the minimum measured volume in the cuvette. The ratio R1: sample does not change.

UNIT CONVERSION

g/dL × 10 = g/L

EXPECTED VALUES⁷

In Serum:

Adults	20–60 years:	3.5–5.2 g/dL	Children	0–4 days	2.8–4.4 g/dL
	60–90 years:	3.2–4.6 g/dL		4 days–14 years	3.8–5.4 g/dL
				14–18 years	3.2–4.5 g/dL

It is recommended that each laboratory verifies this range or derives reference interval for the population it serves.

ANALYTICAL PERFORMANCE

Data contained within this section is representative for performance on ERBA XL-640 automatic system. Data obtained in your laboratory may differ from these values. Data for other ERBA XL automatic systems are available on www.erba.com.

Limit of quantification: 0.07 g/dL

Limit of quantification represents the lowest measurable analyte level. It is calculated as the determined activity of diluted sample to have CV <20 % (n = 30).

Linearity: 6.50 g/dL

Linearity is the highest measured activity with recovery within ±10 % from theoretical value.

Precision:

Precision was determined by using controls in an internal protocol with repeatability (n = 20) and intermediate precision (2 aliquots per run, 2 run per day, 20 days). The following results were obtained:

Repeatability	Mean (g/dL)	SD (g/dL)	CV (%)
Sample 1	3.46	0.016	0.46
Sample 2	4.85	0.025	0.52

Intermediate precision	Mean (g/dL)	SD (g/dL)	CV (%)
Sample 1	3.43	0.079	2.31
Sample 2	4.86	0.066	1.35

Accuracy

Two different validated control materials were used. Determined bias is and -6.1 % at the target value 4.48 g/dL and -8.6 % at the target value 5.26 g/dL.

Comparison

A comparison between XL-640 automatic system ALBUMIN (y) and a commercially available test (x) using 146 samples gave following results:

Linear regression:
 $y = 0.940x + 0.132$ g/dL $r = 0.977$

Passing-Bablok⁸:
 $y = 0.962x + 0.064$ g/dL $r = 0.963$

Interferences

Criterion: Recovery within ±10 % of initial value of albumin concentration in the sample (serum) without interfering substance.
 Following substances do not interfere: haemoglobin up to 4.5 g/L, bilirubin up to 40 mg/dL, triglycerides up to 650 mg/dL.

Limitations:

- Deteriorated reagents (e.g. exceeding the storage temperature) may give incorrect results.
- Quality of reagents is monitored on automatic systems ERBA XL by checking of the maximum permissible absorbance value of blank.
- High concentration of haemoglobin, bilirubin and triglycerides in sample can interfere with determination of albumin. See paragraph Interferences.
- Drugs: No interference was found at therapeutic concentrations using common drug panels⁹.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

Hazards identification in accordance with Regulation (EC) No 1272/2008

R1

UFI: Y8RN-RJU5-0E7S-53CM

Warning

Hazard statement:

H319 Causes serious eye irritation.

Precautionary statement:

P280 Wear protective gloves/protective clothing/eye protection.
 P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
 P337+P313 If eye irritation persists: Get medical advice/attention.

WASTE MANAGEMENT

Please refer to local legal requirements.

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1. Dumas BT, Watson WA, Biggs HG. Albumin standards and the measurement of serum albumin with bromocresol green. Clin Chim Acta 1971;31:87-96.
2. Leonard, P. L., Persaud, J., Motwani, R.: Clin. Chim. Acta 35, 409, 1971.
3. Grant GH, Silverman LM, Christenson RH. Amino acids and proteins. In: Tietz NW, ed. Fundamentals of Clinical Chemistry, 3rd edition Philadelphia, PA: WB Saunders 1987:328-330.
4. Marshall WJ, ed. Illustrated Textbook of Clinical Chemistry, 3rd ed. London: Gower Medical Publishing 1989;207-218.
5. Use of Anticoagulants in Diagnostic Laboratory Investigations. WHO Publication WHO/DIL/LAB/99.1 Rev. 2: Jan 2002.
6. Baudner S, Bienvenu J, Blirup-Jensen S, et al. The certification of a matrix reference material for immunochemical measurement of 14 human serum proteins CRM470. Report EUR 15243 EN 1993;1-186.
7. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. Burtis, C.A., Ashwood, E.R., Bruns, D.E.; 5th edition, WB Saunders Comp., 2012.
8. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11): 783-790.
9. Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 2001;38:376-385.

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