AST/GOT

Cat. No.	Pack Name	Packaging (Content)
XSYS0016	AST/GOT 330	R1: 6 × 44 mL, R2: 6 × 11 mL, RFID tag, instruction for use

INTENDED USE

The kit is intended for *in vitro* photometric quantitative determination of aspartate aminotransferase in human serum and plasma on automatic systems ERBA XL. In combination with other parameters it is intended for screening, monitoring and diagnosis of liver diseases. For professional use in clinical laboratories only.


CLINICAL SIGNIFICANCE

Aspartate aminotransferase (AST/GOT) is widely distributed with high concentrations in the heart, liver, skeletal muscle, kidney and erythrocytes. Damage or disease to any of these tissues such as myocardial infarction, viral hepatitis, liver necrosis, cirrhosis and muscular dystrophy may result in raised levels of AST/GOT. Following myocardial infarction, serum AST is elevated and reaches a peak two days after onset.

In patients undergoing renal dialysis or those with vitamin B6 deficiency, serum AST/GOT may be decreased. The apparent reduction in AST may be related to decreased pyridoxal phosphate, the prosthetic group for AST/GOT, resulting in an increase in the ratio of apoenzyme to holoenzyme. Addition of pyridoxal-5-phosphate, recommended by IFCC¹, stabilizes the activity of AST/GOT and avoids falsely low values in samples containing insufficient endogenous pyridoxal-5-phosphate, e.g. from patients with myocardial infarction, liver disease and intensive care patients.

PRINCIPLE

This assay follows the recommendations of the IFCC¹ without pyridoxal phosphate. AST/GOT in the sample catalyzes the transfer of an amino group between L-aspartate and 2-oxoglutarate to form oxaloacetate and L-glutamate. The oxaloacetate then reacts with NADH, in the presence of malate dehydrogenase (MDH), to form NAD⁺. The series of reactions involved in the assay system is as follows:

The reaction is monitored by measuring the rate of decrease in absorbance at 340 nm due to the oxidation of NADH. The rate of oxidation of NADH is proportional to the activity of AST/GOT present in the sample.

Addition of lactate dehydrogenase to the reagent is necessary to achieve rapid and complete reduction of endogenous pyruvate so that it does not interfere with the assay.

REAGENT DESCRIPTION AND COMPOSITION

R1		R2	
Tris buffer (pH 7.8)	110 mmol/L	CAPSO	20 mmol/L
L-aspartic acid	340 mmol/L	2-oxoglutarate	85 mmol/L
LDH	≥2250 U/L	NADH	1.05 mmol/L
MDH	≥810 U/L	Sodium azide	0.95 g/L
Sodium chloride	200 mmol/L		
Sodium azide	0.95 g/L		

COMPOSITION OF REACTION MIXTURE

Tris buffer	81.5 mmol/L
L-aspartic acid	252 mmol/L
LDH	≥1667 U/L
MDH	≥600 U/L
CAPSO	3.7 mmol/L
2-oxoglutarate	15.7 mmol/L
NADH	0.19 mmol/L
Sodium chloride	148 mmol/L
Sodium azide	0.88 g/L

REAGENT PREPARATION

Reagents are liquid, ready to use. Load the number of tests from the RFID tag before using a new kit.

MATERIAL REQUIRED BUT NOT PROVIDED WITH THE DEVICE

XL MULTICAL 4×3, Cat. No. XSYS0034
 XL MULTICAL 10×3, Cat. No. XSYS0122
 ERBA NORM 4×5, Cat. No. BLT00080
 ERBA NORM 10×5, Cat. No. XSYS0123
 ERBA PATH 4×5, Cat. No. BLT00081
 ERBA PATH 10×5, Cat. No. XSYS0124
 Erba XL analysers: XL-200, Cat. No. INS00002
 XL-640, Cat. No. INS00008
 XL-1000, Cat. No. INS00010

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C. On board stability: min. 30 days if refrigerated 2–10 °C and not contaminated.

SPECIMEN COLLECTION AND HANDLING

It is recommended to follow ISO 15189 and laboratory instruction.

For specimen collection and preparation only use suitable tubes or collection containers.

Only the specimens listed below were tested and found acceptable.

Serum.

Plasma: Li-heparin and K₂-EDTA plasma.

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer. Centrifuge samples containing precipitates before performing the assay.

See the limitations and interferences section for details about possible sample interferences.

Stability in serum / plasma⁴:	4 days at	15–25 °C
	7 days at	2–8 °C
	3 months at	-20 °C

Discard contaminated specimens.

CALIBRATION

Calibration with calibrator XL MULTICAL is recommended.

2 point calibration (blank and calibrator); distilled water is recommended as blank

Calibration frequency: 30 days

Calibration is needed:

- after reagent lot change
- as required by internal quality control procedures
- calibration interval may be extended based on acceptable verification of calibration by the laboratory.

QUALITY CONTROL

For quality control ERBA NORM and ERBA PATH are recommended. The control intervals and limits should be adapted according to each individual laboratory's requirements. Values obtained should fall within the defined intervals. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

TRACEABILITY

This method, calibrator XL MULTICAL and controls ERBA NORM and ERBA PATH have been standardized to IFCC recommended method¹.

ASSAY PROCEDURE AND CALCULATION

ERBA XL automatic systems calculate the concentration of each sample. For assay parameters see www.erba.com.

Assay parameters for ERBA XL automatic systems

Assay type	Rate A
Curve type	Linear
Wavelength (prim. / sec.)	340 / 660 nm
Reading time	60–240 s after adding of R2
Reaction direction	Decrease
Unit	U/L (μkat/L)
Reagent volumes	
R1	160 μL
R2	40 μL
Sample	16 μL

Note: reagents and sample volumes can be different for individual ERBA XL automatic systems depending on the minimum measured volume in the cuvette. The ratio R1:R2:sample does not change.

UNIT CONVERSION

U/L × 0.0167 = μkat/L

EXPECTED VALUES⁵

At 37 °C

Females <31 U/L

Males <35 U/L

It is recommended that each laboratory verifies this range or derives reference interval for the population it serves.

ANALYTICAL PERFORMANCE

Data contained within this section is representative for performance on ERBA XL-640 automatic system. Data obtained in your laboratory may differ from these values. Data for other ERBA XL automatic systems are available on www.erba.com.

Limit of quantification: 2.13 U/L

Limit of quantification represents the lowest measurable analyte level. It is calculated as the determined activity of diluted sample to have CV <20 % (n = 30).

Linearity: 418 U/L

Linearity is the highest measured activity with recovery within ±10 % from theoretical value.

Precision:

Precision was determined by using controls in an internal protocol with repeatability (n = 20) and intermediate precision (2 aliquots per run, 2 run per day, 20 days). The following results were obtained:

Repeatability	Mean (U/L)	SD (U/L)	CV (%)
Sample 1	42.4	0.60	1.41
Sample 2	153.3	1.13	0.74

Intermediate precision	Mean (U/L)	SD (U/L)	CV (%)
Sample 1	42.5	1.15	2.71
Sample 2	151.7	2.06	1.36

Accuracy

Two different validated control materials were used. Determined bias is -4.8 % at the target value 146.6 U/L and -4.2 % at the target value 197.2 U/L.

Comparison

A comparison between XL-640 automatic system AST/GOT with PDP (y) and a commercially available test (x) using 48 samples gave following results:

Linear regression:

y = 1.025x - 2.74 U/L r = 0.999

Passing-Bablok⁶:

y = 1.000x - 1.80 U/L r = 0.977

Interferences

Criterion: Recovery within ±10 % of initial value of AST/GOT activity in the sample without interfering substance.

Following substances do not interfere: bilirubin up to 40 mg/dL, triglycerides up to 850 mg/dL. Haemolysis interferes due to AST/GOT activity from erythrocytes.

Drugs: No interference was found at therapeutic concentrations using common drug panels⁷.

Limitations:

- Deteriorated reagents (e.g. exceeding the storage temperature) may give incorrect results. Quality of reagents is monitored on automatic systems ERBA XL by checking of the minimal permissible absorbance value of blank.

- High concentration of haemoglobin, bilirubin and triglycerides in sample can interfere with determination of AST/GOT. See paragraph Interferences.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

Hazards identification in accordance with Regulation (EC) No 1272/2008

R1

UFI: YMDF-7EN4-1H7J-XF4E

Warning

Contains: aspartic acid, tris(hydroxymethyl)aminomethane, sodium hydroxide

Hazard statement:

H315 Causes skin irritation.

H319 Causes serious eye irritation.

Precautionary statement:

P280 Wear protective gloves / protective clothing / eye protection.

P303+P352 IF ON SKIN: Wash with plenty of water and soap.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P337+P313 If eye irritation persists: Get medical advice/attention.

R2

The reagent is not classified as dangerous.

WASTE MANAGEMENT

Please refer to local legal requirements.

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1. Schumann G, Bonora R, Ceriotti F, Féraud G et al. IFCC primary reference procedure for the measurement of catalytic activity concentrations of enzymes at 37 °C. Part 5: Reference procedure for the measurement of catalytic concentration of aspartate aminotransferase. Clin Chem Lab Med 2002;40:725-33.
2. Thomas L. Alanine aminotransferase (ALT), Aspartate aminotransferase (AST). In: Thomas L, editor, Clinical Laboratory Diagnostics, 1st ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998, p. 55-65.
3. Bergmeyer HU, Horder M, Rej R, Approved Recommendation (1985) on IFCC Methods for the Measurement of Catalytic Concentration of Enzymes, L. Clin. Chem. Clin. Biochem. 1986; 24: 481-495
4. Use of Anticoagulants in Diagnostic Laboratory Investigations. WHO Publication WHO/DIL/LAB/99.1 Rev. 2: Jan 2002.
5. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. Burtis, C.A., Ashwood, E.R., Bruns, D.E.; 5th edition, WB Saunders Comp., 2012.
6. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11):783-790.
7. Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 2001;38:376-385

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