
CREATININE

Cat. No.	Pack Name	Packaging (Content)
XSYS0024	CREA 275	R1: 5 × 44 mL, R2: 5 × 11 mL, RFID tag, instruction for use, chimneys: 2 pcs

INTENDED USE

The kit is intended for *in vitro* photometric quantitative determination of creatinine in human serum, plasma and urine on automatic systems ERBA XL. In combination with other parameters it is intended for screening, monitoring of renal function, diagnosis of renal diseases. For professional use in clinical laboratory only.

CLINICAL SIGNIFICANCE

Creatinine is a waste product excreted by the kidneys mainly by glomerular filtration. The concentration of creatinine in plasma of a healthy individual is fairly constant, independent from water intake, exercise and rate of urine production. Therefore, increased plasma creatinine values always indicate decreased excretion, i.e. impaired kidney function. Creatinine clearance is a good indicator for the glomerular filtration rate (GFR) which allows better detection of kidney diseases and monitoring of renal function. For this purpose, creatinine is measured simultaneously in serum and urine collected over a defined time period. The serum creatinine levels do not start to rise until renal function has decreased by at least 50 %.

PRINCIPLE

Creatinine reacts with alkaline picrate to produce an orange-red complex (Jaffe reaction). This is a non-specific reaction and is given by many other substances. Specificity of the assay has been improved by the introduction of a kinetic method, however, the cephalosporin antibiotics are still major interferents^{1,2,3,4,5,6}.

The rate of complex formation, measured at 505 nm, is proportional to the creatinine concentration in the sample.

REAGENT DESCRIPTION AND COMPOSITION

R1		R2	
Sodium hydroxide	240 mmol/L	Picric acid	26 mmol/L
COMPOSITION OF REACTION MIXTURE			
Sodium hydroxide	183 mmol/L		
Picric acid	5 mmol/L		

REAGENT PREPARATION

Reagents are liquid, ready to use. Load the number of tests from the RFID tag before using a new kit.

MATERIAL REQUIRED BUT NOT PROVIDED WITH THE DEVICE

XL MULTICAL 4×3, Cat. No. XSYS0034
XL MULTICAL 10×3, Cat. No. XSYS0122
ERBA NORM 4×5, Cat. No. BLT00080
ERBA NORM 10×5, Cat. No. XSYS0123
ERBA PATH 4×5, Cat. No. BLT00081
ERBA PATH 10×5, Cat. No. XSYS0124
Erba XL analysers: XL-200, Cat. No. INS00002
XL-640, Cat. No. INS00008
XL-1000, Cat. No. INS00010

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.

On board stability: min. 30 days if refrigerated (2–10 °C) and not contaminated and if chimneys are used. Absorption of atmospheric CO₂ in open reagent bottles R1 and R2 leads to the reduced stability of reagents. To minimize the rate of absorption of CO₂, the use of "chimneys" provided with the kit is recommended. Use of chimneys has shown improvement in calibration stability, provided the storage and use of reagents is as per instructions.

Chimneys provided with the kit are to be placed directly into the reagent bottle immediately after opening the cap. White chimney is to be used for reagent R1 and black chimney for reagent R2. Colour differentiation of chimneys is just to ensure chimneys are not exchanged between R1 & R2. Chimneys can be re-used again for the bottle of the same kit. Before using the chimney in a new bottle, it is recommended to wash chimneys well with distilled water and dry them thoroughly.

SPECIMEN COLLECTION AND HANDLING

It is recommended to follow ISO 15189 and laboratory instruction.

For specimen collection and preparation only use suitable tubes or collection containers.

Only the specimens listed below were tested and found acceptable.

Serum.

Plasma: Li-heparin and K₂-EDTA plasma.

Urine: Collect urine without using additives. If urine must be collected with a preservative for other analytes, only hydrochloric acid (14 to 47 mmol/L urine, e.g. 5 mL 10 % HCl or 5 mL 30 % HCl per liter urine) or boric acid (81 mmol/L, e.g. 5 g per liter urine) may be used. Dilute urine samples using redistilled water in 1+19 ration and multiply results by 20.

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay.

See the Limitations and Interferences section for details about possible sample interferences.

Stability in serum / plasma ² :	
7 days at 15–25 °C	
7 days at 2–8 °C	
3 months at -20 °C	
Stability in urine ² :	
2 days at 15–25 °C	
6 days at 2–8 °C	
6 months at -20 °C	

Discard contaminated specimens.

CALIBRATION

Calibration with calibrator XL MULTICAL is recommended.

2 point calibration (blank and calibrator); distilled water is recommended as blank

Calibration frequency: 7 days

Calibration is needed:

- after reagent lot change
- as required by internal quality control procedures
- calibration interval may be extended based on acceptable verification of calibration by the laboratory

QUALITY CONTROL

For quality control ERBA NORM and ERBA PATH are recommended.

The control intervals and limits should be adapted according to each individual laboratory's requirements. Values obtained should fall within the defined intervals. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

TRACEABILITY

This method, calibrator XL MULTICAL and controls ERBA NORM and ERBA PATH have been standardized against to ID/MS.

ASSAY PROCEDURE AND CALCULATION

ERBA XL automatic systems calculate the concentration of each sample. For assay parameters see www.erba.com.

Assay parameters for ERBA XL automatic systems

Assay type	Rate A
Curve type	Linear
Wavelength (prim./sec.)	505/570 (578) nm
Reading time	30–150 s after adding of R2
Reaction direction	Increasing
Unit	mg/dL (μmol/L)
Reagent volumes	
R1	160 μL
R2	40 μL
Sample volumes	10 μL

Note: reagents and sample volumes can be different for individual ERBA XL automatic systems depending on the minimum measured volume in the cuvette. The ratio R1:R2:sample does not change.

NOTE

The determination of serum/plasma creatinine passed on Jaffe method is not fully specific. Alkaline picrate reacts also with some others substances present in serum matrix. Therefore, to correct the matrix effects, it is recommended to use correction factor in ASSAY PARAMETERS:

$$y = ax + b$$

$$b = -0.2 \text{ mg/dL} / b = -18 \text{ μmol/L}$$

UNIT CONVERSION

$$\text{mg/dL} \times 88.4 = \text{μmol/L}$$

EXPECTED VALUES

Serum¹⁰

Cord	0.6–1.2 mg/dL	18–60 y: Male	0.9–1.3 mg/dL
Newborn (1–4 d)	0.3–1.0 mg/dL	Female	0.6–1.1 mg/dL
Infant	0.2–0.4 mg/dL	60–90 y: Male	0.8–1.3 mg/dL
Child	0.3–0.7 mg/dL	Female	0.6–1.2 mg/dL
Adolescent	0.5–1.0 mg/dL	>90 y: Male	1.0–1.7 mg/dL
		Female	0.6–1.3 mg/dL

Urine¹⁰

Infant	8–20 mg/kg/day	Adult: Male	14–26 mg/kg/day
Child	8–22 mg/kg/day	Female	11–20 mg/kg/day
Adolescent	8–30 mg/kg/day		

It is recommended that each laboratory verifies this range or derives reference interval for the population it serves.

ANALYTICAL PERFORMANCE

Data contained within this section is representative for performance on ERBA XL-640 automatic system. Data obtained in your laboratory may differ from these values. Data for other ERBA XL automatic systems are available on www.erba.com.

Limit of quantification:

Serum/plasma	0.020 mg/dL
Urine	0.38 mg/dL

Limit of quantification represents the lowest measurable analyte level. It is calculated as the determined activity of diluted sample to have CV <20 % (n = 30).

Linearity:

Serum/plasma	30.0 mg/dL
Urine	600 mg/dL

Linearity is the highest measured activity with recovery within ±10 % from theoretical value.

Precision:

Precision was determined by using controls in an internal protocol with repeatability (n = 20) and intermediate precision (2 aliquots per run, 2 run per day, 20 days). The following results were obtained:

Repeatability (serum)	Mean (mg/dL)	SD (mg/dL)	CV (%)	Repeatability (urine)	Mean (mg/dL)	SD (mg/dL)	CV (%)
Sample 1	1.34	0.015	1.10	Sample 1	56.2	0.73	1.31
Sample 2	3.34	0.042	1.24	Sample 2	125.2	1.70	1.36

Intermediate precision (serum)	Mean (mg/dL)	SD (mg/dL)	CV (%)	Intermediate precision (urine)	Mean (mg/dL)	SD (mg/dL)	CV (%)
Sample 1	1.40	0.042	3.02	Sample 1	61.8	1.89	3.05
Sample 2	3.56	0.129	3.62	Sample 2	134.3	4.07	3.03

Accuracy

Two different validated control materials for serum and urine were used. Determined bias is -1.9 % at the target value 1.38 mg/dL, -7.5 % at the target value 3.99 mg/dL for serum, -12.2 % at the target value 66.9 mg/dL and -14.2 % at the target value 141.4 mg/dL for urine.

Comparison

A comparison between XL-640 automatic system CREATININE (y) and a commercially available test (x) using 149 samples (serum) gave following results:

Linear regression:

$$y = 0.955x + 0.131 \text{ mg/dL} \quad r = 0.996$$

Passing-Bablok¹²:

$$y = 0.973x + 0.111 \text{ mg/dL} \quad r = 0.992$$

Interferences

Criterion: Recovery within ±10 % of initial value of creatinine concentration in the sample (serum) without interfering substance.

Following substances do not interfere: haemoglobin up to 10 g/L, bilirubin up to 4 mg/dL, triglycerides up to 850 mg/dL.

Cefoxitin, Cephalothin and Cyanokit (Hydroxocobalamin) may cause false-positive results in serum or plasma^{13,14,15,16}.

Limitations:

- Deteriorated reagents (e.g. exceeding the storage temperature) may give incorrect results. Quality of reagents is monitored on automatic systems ERBA XL by checking of the maximal permissible absorbance value of blank.
- High concentration of haemoglobin, bilirubin and triglycerides in sample can interfere with determination of creatinine. See paragraph Interferences.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

Hazards identification in accordance with Regulation (EC) No 1272/2008

R1

UFI: KWFQ-3J78-FE72-MFJX

Warning

Contains: sodium hydroxide

Hazard statement:

H315 Causes skin irritation.

H319 Causes serious eye irritation.

Precautionary statement:

P280 Wear protective gloves / protective clothing / eye protection.

P302 + P352 IF ON SKIN: Wash with plenty of water.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

R2

Reagent is not classified as dangerous.

WASTE MANAGEMENT

Please refer to local legal requirements.

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1. Jaffé M. Ueber den Niederschlag, welchen Pikrinsäure in normalem Harn erzeugt und über eine neue Reaktion des Kreatinins. Z Physiol Chem 10: 391-400, 1886.
2. Fabiny DL, Ertinghausen G. Automated reaction-rate method for determination of serum creatinine with the CentrifChem Clin Chem. 17: 696-700, 1971.
3. Bartels H, Böhmer M. Micro-determination of creatinine. Clin Chim Acta 32: 81-85, 1971.
4. Myers, G. L., Greg Miller, W., Coresh, J., Fleming, J., Greenberg, N. et al.: Recommendations for Improving Serum Creatinine Measurement, Clin. Chem. 52: 5-18, 2006.
5. Fridecký B. Program zlepšování kvality měření sérového kreatininu, Klin. Biochem. Metab., 14 (35), No.3: 173-176, 2006.
6. Bowers, L. D., Wong, E. T. Kinetic serum creatinine assays. II. A critical evaluation and review, Clin. Chem. 26: 555, 1980.
7. Guder W, Fonseca-Wollheim W, Ehret W, et al. Die Qualität Diagnostischer Proben, 6. Aufl. Heidelberg: BD Diagnostics, 2009.
8. Newman DJ, Price CP. Renal function and nitrogen metabolites. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company: 1204-1270, 1999.
9. Thomas L. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft: 366-374, 1998.
10. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. Burtis CA. Ashwood. ER. Bruns. DE; 5th edition. WB Saunders Comp.. 2012.
11. Tietz Fundamentals of Clinical Chemistry and Molecular Diagnostics. Rifai N. Horvath AR. Wittwer CT; 8th edition. Elsevier. 2019.
12. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem, Nov;26(11): 783-790, 1988.
13. Sonntag O, Scholer A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 38, 376-385, 2001.
14. Kroll MH, Elin RJ, Mechanism of cefoxitin and cephalothin interference with the Jaffé method for creatinine, Clin Chem 29, 12, 2044-2048, 1983.
15. Wang P, Commentary on Impact of Cyanokit Administration on Point-of-Care Glucose and Laboratory Test Accuracy, Clin Chem 68, 4, 509-510, 2022.
16. Green AJE, Halloran SP, Mould GP, Barbour HM, Pritchard JL, Halloran MJ, Labib M, Interference by newer cephalosporins in current methods for measuring creatinine, Clin Chem 36, 12, 2139-2140, 1990.

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Dátum revízie: 16. 12. 2025

LITERATÚRA

1. Jaffé M. Ueber den Niederschlag, welchen Pikrinsäure in normalem Harn erzeugt und über eine neue Reaktion des Kreatinins. Z Physiol Chem 10: 391-400, 1886.
2. Fabiny DL, Ertinghausen G. Automated reaction-rate method for determination of serum creatinine with the CentrifChem Clin Chem. 17: 696-700, 1971.
3. Bartels H, Böhmer M. Micro-determination of creatinine. Clin Chim Acta 32: 81-85, 1971.
4. Myers, G. L., Greg Miller, W., Coresh, J., Fleming, J., Greenberg, N. et al.: Recommendations for Improving Serum Creatinine Measurement, Clin. Chem. 52: 5-18, 2006.
5. Fridecký B. Program zlepšování kvality měření sérového kreatininu, Klin. Biochem. Metab., 14 (35), No.3: 173-176, 2006.
6. Bowers, L. D., Wong, E. T. Kinetic serum creatinine assays. II. A critical evaluation and review, Clin. Chem. 26: 555, 1980.
7. Guder W, Fonseca-Wollheim W, Ehret W, et al. Die Qualität Diagnostischer Proben, 6. Aufl. Heidelberg: BD Diagnostics, 2009.
8. Newman DJ, Price CP. Renal function and nitrogen metabolites. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company: 1204-1270, 1999.
9. Thomas L. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft: 366-374, 1998.
10. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. Burtis CA. Ashwood. ER. Bruns. DE; 5th edition. WB Saunders Comp.. 2012.
11. Tietz Fundamentals of Clinical Chemistry and Molecular Diagnostics. Rifai N. Horvath AR. Wittwer CT; 8th edition. Elsevier. 2019.
12. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem, Nov;26(11): 783-790, 1988.
13. Sonntag O, Scholter A. Drug interference in clinical chemistry: recommendation of drugs and their concentrations to be used in drug interference studies. Ann Clin Biochem 38, 376-385, 2001.
14. Kroll MH, Elin RJ, Mechanism of cefoxitin and cephalothin interference with the Jaffé method for creatinine, Clin Chem 29, 12, 2044-2048, 1983.
15. Wang P, Commentary on Impact of Cyanokit Administration on Point-of-Care Glucose and Laboratory Test Accuracy, Clin Chem 68, 4, 509-510, 2022.
16. Green AJE, Halloran SP, Mould GP, Barbour HM, Pritchard JL, Halloran MJ, Labib M, Interference by newer cephalosporins in current methods for measuring creatinine, Clin Chem 36, 12, 2139-2140, 1990.

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Katalógové číslo

Číslo šarže

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eIFU:
www.erba.com

Diagnostický zdravotnícký prostriedok *in vitro*


Výrobca

Obmedzenie teploty

Obsah

Erba Lachema s.r.o., Karásek 2219/1d, 621 00 Brno, CZ
e-mail: diagnostics@erba.com, www.erba.com

CC/IFU/058/25/A

Dátum revízie: 16. 12. 2025