HDL/LDL CAL

Cat. No.	Pack Name	Packaging (Content)
XSYS0061	HDL/LDL CAL	2 × 1 mL, instruction for use

INTENDED USE

HDL/LDL CAL is calibration material for *in vitro* photometric quantitative determination of HDL/LDL cholesterol in human serum and plasma on automatic systems ERBA XL. Intended for screening, monitoring and diagnosis of coronary heart disease and lipoprotein metabolism disorders. For professional use in clinical laboratory only.

PRINCIPLE

HDL/LDL CAL is a lyophilized calibrator based on human serum. The concentrations of individual analytes are suitable for calibrations of clinical chemistry assays on automated clinical chemistry analysers Erba XL.

DESCRIPTION AND COMPOSITION

The calibrator is based on human serum with chemical additives and stabilisers; it is provided in lyophilized form for increase stability.

CERTIFIED VALUE

Calibrator LOT	Expiry	Parameter	Concentration	
xxxxxxx	yyyy-mm	HDL	00.0 mg/dL	0.00 mmol/L
		LDL	00.0 mg/dL	0.00 mmol/L

UNIT CONVERSION

mg/dL × 0.026 = mmol/L

RECONSTITUTION

Open the bottle very carefully and add exactly 1.00 ml of redistilled or deionized water. Reconstitute the calibrator within 30 min at room temperature of (+15 to +25) °C in a dark place. Swirl the bottle carefully, do not shake (avoid foaming)!

STORAGE AND HANDLING

Unreconstituted lyophilized HDL/LDL CAL is stable till the expiry date mentioned on the label when stored at 2–8 °C.

The reconstituted calibrator is stable for at least:

7 days at (+2 to +8) °C in dark

Data concerning calibrator stability are valid only for the calibrator not contaminated microbially. Do not use the calibrator with a turbidity developed pointing at a possible bacterial growth.

MATERIALS REQUIRED, BUT NOT PROVIDED

HDL C 160, Cat. No. XSYS0043
HDL 80, Cat. No. BLT00028
LDL C 80, Cat. No. XSYS0044
LDL 80, Cat. No. BLT00041
ERBA NORM 4×5, Cat. No. BLT00080
ERBA NORM 10×5, Cat. No. XSYS0123
ERBA PATH 4×5, Cat. No. BLT00081
ERBA PATH 10×5, Cat. No. XSYS0124
Erba XL analysers: XL-200, Cat. No. INS00002
XL-640, Cat. No. INS00008
XL-1000, Cat. No. INS00010

or any instrument with temperature control of 37 ±0.5 °C that is capable of reading absorbance at 600/700 nm

ASSAY PROCEDURE

Use calibrator as specified in the relevant Instruction For Use for the Erba reagents.

EVALUATION

The XL-analysers perform the evaluation automatically.

EXPECTED RESULTS

The calibrator values were determined using the method stated in the enclosed value sheets. Determinations were performed under strictly standardized conditions on XL analysers using Erba system reagents and the standard reference materials.

The calibrator values were obtained via single determinations performed in different laboratories in several independent series. The calibrator value specified is the mean of all values obtained.

Calibrator was tested within the functional testing of all diagnostic kits.

TRACEABILITY

This calibrator has been standardized against the reference material NIST SRM 1951.

PERFORMANCE CHARACTERISTICS

Performance characteristics of individual analytes on the Erba XL automatic systems are available on websites www.erba.com.

LIMITATIONS

Erroneous results can occur from reconstitution inaccuracies and technical errors associated with the assay procedure.

Do not use the reconstituted calibrator if there is visible evidence of microbial growth in the vial.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

Health protection

Serum used in the manufacture of this product has been tested and found non-reactive for hepatitis B surface antigen (HbsAg), antibody to Hepatitis C (HCV), HIV-1 p24 antigen, and antibody to HIV-1/2. The test procedures do not guarantee that all infectious agents will be detected. Because no test method can offer complete assurance that Hepatitis B virus Hepatitis C virus and HIV-1/2 or other infectious agents are absent, the material should be handled as potentially infectious. During the work all principles of general safety as well as the principles of work with biological materials are to be maintained. Above all observe the principles of personal hygiene, use personal protection aids (clothes, shoes, etc.), do not eat, drink and smoke. After finishing the work and before a meal wash your hands with soap and warm water and apply a suitable regenerating cream. When the bottle is broken, the contaminated surface must be disinfected with Persteril or Chloramine.

Hazards identification in accordance with Regulation (EC) No 1272/2008

The calibrator is not classified as dangerous.

WASTE MANAGEMENT

Please refer to local legal requirements.

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1. *In vitro* diagnostic medical devices – Requirements for establishing metrological traceability of values assigned to calibrators, trueness control materials and human samples. ISO 17511, 2nd edition, 2020.
2. *In vitro* diagnostic medical devices – Requirements for international harmonisation protocols establishing metrological traceability of values assigned to calibrators and human samples. ISO 21151, 1st edition, 2021.

USED SYMBOLS / POUŽITÉ SYMBOLY / УСЛОВНЫЕ ОБОЗНАЧЕНИЯ / POUŽITÉ SYMBOLY / SÍMBOLOS UTILIZADOS ВИКОРИСТАНІ ПОЗНАЧКИ / SYMBOLES UTILISÉS / SÍMBOLOS USADOS

REF	Catalogue number Katalogové číslo Номер по каталогу Katalógové číslo Número de catálogo Каталожний номер Numéro de catalogue Número de catálogo	LOT	Lot number Číslo šarže Код партии Číslo šarže Número de lote Номер партії Numéro de lot Número de lote		Expiry date Datum expirace Использовать до Dátum expirácie Fecha de caducidad Термін придатності Date d'expiration Data de validade		Consult instructions for use Čtěte návod k použití Обратитесь к инструкции по применению или к инструкции по применению в электронном виде Čítajte návod k použitiu Consulte las instrucciones de uso Перед використанням уважно вивчіть інструкцію Consulter la notice d'utilisation Veja as instruções de uso
IVD	<i>In vitro</i> diagnostic medical device Diagnostický zdravotnický prostředek <i>in vitro</i> Медицинское изделие для диагностики <i>in vitro</i> Diagnostický zdravotnícky prostriedok <i>in vitro</i> Dispositivo médico para diagnóstico <i>in vitro</i> <i>In vitro</i> діагностика Dispositif médical de diagnostic <i>in vitro</i> Diagnóstico <i>in vitro</i>		Manufacturer Výrobce Изготовитель Výrobca Fabricante Виробник Fabricant Fabricante		Temperature limit Omezení teploty Температурный диапазон Obmedzenie teploty Límite de temperatura Температура зберігання Limites de température Temperatura de armazenamento	CONT	Content Obsah Содержание Obsah Contenido Вміст Contenu Conteúdo
	Biological risks Biologická rizika Биологические риски Biologické riziká Riesgos biológicos Біологічні ризики Risques biologiques Riscos biológicos						